



New Concepts, New Connections™



1-866-305-6145

How to Fit




Putting On

1. Open snap buckle.
2. Place RIGHT hand over the harness label on B strap.
3. Face your dog and place B strap over his head and onto his back.
The label will be on your dog's left side and the C strap will be over your dog's front (forechest).

Fitting

Adjust (4) straps in the following order:

- 1 B Strap
- 2 G Strap
- 3 C Strap (both sides)



1. B Strap (important)

- Adjust the B strap to raise side rings.
- Place side rings in the middle of your dog's sides (between your dog's back and girth).

2. G Strap

- Adjust the G strap to a snug fit (like a comfortable belt).
- Reach under your dog with your LEFT hand and grasp the G strap to buckle.
- Check for comfort room (one finger under strap) when your dog sits. The G strap needs to be snug to prevent the side rings from shifting forward (over your dog's shoulders) but not tight so your dog can sit comfortably.

3. C Strap


- Adjust both sides of C strap (carefully) to place the leash ring over your dog's breast bone.
- **Do not over-adjust.** It is better to have the C strap loose than too tight.
- Have your dog sit. If the C strap is horizontal with the leash ring held up by your dog's breast bone, you have adjusted the C strap enough.
- The C strap will be horizontal **only** when your dog sits or when held up by your leash. When you look at your dog from the side while he stands, you should see the side rings are higher than the leash ring.
- Your dog should mainly feel the B & G straps (belt behind front legs).

Check The Fit

- The Side Rings should be positioned in the middle of your dog's sides.
- The Chest Strap is **not** horizontal when your dog stands (side rings are higher than the leash ring).
- Your dog's chest expands when he sits and naturally raises the chest strap to a horizontal position.



Correct Fit



Wrong Fit

Wrong Fit

- The C strap is **too low and too tight**. This restricts the dog's shoulders and front legs.
- The B strap should be shortened, the G strap lengthened, and the C strap(s) lengthened.

[Download our Fitting Guide](#)




New Concepts, New Connections™



1-866-305-6145

Measurement Instruction



Girth Measurement

Using a flexible measuring tape (or leash), measure from just behind the front legs all the way around your dog. You are measuring the circumference of the body just behind the front legs, like measuring for a "belt." Then go to the Products page (SENSE-ible® or SENSE-ation®) and select the proper size using the girth measurement.

Please note that the measurement is the best way to get the correct size - DO NOT USE THE WEIGHT RANGE to select size.

A = Measuring tape position (just behind front legs)

B = Start and finish at this point


Chest Measurement

The chest measurement is typically 1/2 the girth measurement and is not required for most dogs.

Some breeds, such as Bulldogs, Bull Terriers, etc., will require a measurement as per the diagram (right) to determine a proper fit.

Also, this measurement is necessary for any alteration or modification requests.

Measure around the center of the forechest, over the breastbone. Start behind one of the dog's front legs to the same point behind the other front leg.



Copyright © 2003 - 2015 Softouch Concepts, Inc. All Rights Reserved.
U.S. Patents: #7,107,939 #7,387,087 #7,886,699 #8,869,750.